

MINIMUM THICKNESS




Malleo Dynastab® can be worn in all types of shoes

MAXIMUM IMMOBILISATION



- Resistant and comfortable light metal reinforcement
- Strapping straps

OPTIMUM COMFORT



- Anatomical design
- Micro-ventilated 3D fabric

MALLEO DYNASTAB



INDICATIONS

Post-Trauma

- sprains of the lateral or medial tendons of the ankle due to in- or eversion.

Preventative

- High impact sports that cause high of tension of movement & lateral force on the ankle joint.
eg. squash, netball , rugby.

Post-operative

- Post-moonboot (fracture or partial fractures) to provide stability of the ankle due to weakened tendons

Drop Foot

FEATURES & BENEFITS

Reinforcing wire is designed to fit around the malleoli, (not sub-talar).

- Provides strong stability, yet allows for slight in- and eversion.
- Comfortable - no sub-talar irritation or heel pain.
- Great range of motion possible for plantar and dorsal flexion.

Removable Wire frame for:

- Customising eg bend to accommodate for swollen malleoli.
- Washing purposes

Open heel & Achilles Tendon, preventing chafing of skin on the Achilles.

Stretchable straps with Velcro attachments

- Attached above the Achilles – promotes flexed resting position (aids drop foot management)
- Figure-of-8 binding for little to strong stability.

Laces allow for "easy-on-easy-off" fitting.

Low profile to fit comfortably in almost any shoe or sports shoe.

Material is fully washable and anti-microbial.

Universal

3 Sizes